Zirconia Burs



ジルコニアバー

The grinding of the resin material.

The awesome white blade of zirconia



- Hybrid resin
- Composite resin
- Tech

- Denture Base Resin
- Thermoplastic resin
- Stent



Slippery finished surface

Blade back (the surface on which the resin shavings touches the bur after grinding) with a carefully polished one by one. Smoothly discharge the resin shavings and heat.

"Technology"

Zirconia has good slippery, difficult to hold the frictional heat.

Grinding with high efficiency, and the surface is beautiful.

"Technology"

Beautiful cutting surface

We taking advantage of the many years of carbide bur manufacturing know-how, and the pursuit of grinding efficiency and durability, has developed an original blade shape for zirconia.


Different resin shavings

Compared to resin shabings was grinding with a metal bur, the shavings are not charged with static electricity.

Therefore, it can be removed easily shavings only air from the clothes.



White with cleanliness.





Almost

DD6-Zr20 / DD6-Zr28

Impressive white

The blade made of zirconia. Without causing a blade to be baked with shavings, it can easily be cleaned with a brush or the like.

Zirconia Burs

Part No. : DD6-Zr20 (Fine)

DD6-Zr28 (Super Fine)

Packaging and price:

• 1 pcs = \pm 6,500

(excluding tax)

JAN code:

DD6-Zr20 1 pcs = 4571243041304 DD6-Zr28 1 pcs = 4571243041311

550 5.25 5 500 1572 150 1252

Before use, please read the package insert of the product.



7 – 5 Tsuto-Otsuka Nishinomiya 663 – 8241 Hyogo, JAPAN Tel: +81-798–37–0270 Fax: +81-798–37–0230 www. mt-chuo .jp